


# CODAN Check Valves



The versatile program  
of check valves





# CODAN Check Valves

## **Integrated safety**

CODAN Check Valves are a small detail with maximum importance for patient safety. As a firm or optional integral part in various CODAN products, the check valves reliably prevent reflux or back flow of infusion solution without additional measures and also ensure that blood does not flow back from the patient into the infusion system.

## **Small, perfect and safe**


### **Top-class engineering in a mini-format**

As with all CODAN products, the CODAN Check Valves also manifest decades of experience in research and development.

Precision in production ensures maximum reliability and ease of movement of the CODAN Check Valve. The valve opens upon the slightest conveying pressure. It closes automatically as soon as a counter-pressure results or the conveying pressure has dropped again. Even highly viscous infusions pass through the valves securely and reliably.

The CODAN Check Valves have a particularly low priming volume of 0,15 ml and can be integrated into all infusion systems without any problems.

Be it in standard infusion therapy on wards, in the operating theatre, in intensive care or in anaesthetics – CODAN Check Valves can be used in all areas of everyday hospital work.



# Versatility with precision

## For use in many versions


The matching solution can be found for each application. The CODAN Check Valve can be integrated into an existing system as an individual component. In addition, a variety of connection sets, extensions, connectors, completely pre-fabricated systems and special systems with CODAN Check Valves are available. All variants are available with differing tubing diameters. Connectors with CODAN Check Valves are available as double, triple or quadruple connectors.

## Connections


CODAN Check Valves in CODAN infusion set connections make handling particularly safe and simple – the less individual components are used, the more quickly and more reliably the infusion can be prepared and administered. Special importance is given to the Y-shaped injection piece in safe and low-risk administration of cytostatics.

## Connection variants


The connections provided with a rotating collar which enable secure closure of the Male Luer-Lock without tube torsion are particularly comfortable and user-friendly. Of course CODAN also provides all standard infusion sets with integrated CODAN Check Valves.




R-LOCK,  
Female/Male  
Luer-Lock adapter  
with CODAN  
Check Valve




Inline adapter  
with CODAN  
Check Valve




Double Y-site  
with CODAN  
Check Valves




Y-injection site with  
CODAN Check Valve



Double connector  
with CODAN  
Check Valves



Rotating  
Male Luer-Lock adapter  
with CODAN Check Valve



Male Luer-Lock adapter  
with CODAN Check Valve



CODAN

# CODAN Worldwide



CODAN is known internationally as a manufacturer and supplier of disposable medical transfer systems. The CODAN Companies have more than 1500 employees around the world.

The name CODAN is synonymous with reliability, quality and precision based on the know-how and experience gained from more than 60 years of research and development. Company-owned production facilities and sales companies around the world are a guarantee for efficient production, a tight-knit sales network and a first-class service for customers in the healthcare sector.

## CODAN Companies

CODAN Medizinische Geräte GmbH & Co KG · Deutschland  
CODAN pvb Critical Care GmbH · Deutschland  
CODAN pvb Medical GmbH · Deutschland  
CODAN PORTUGAL, S.A. · Portugal  
CODAN 11, S.A. · Portugal  
CODAN US Corporation · California · USA  
CODAN Inc. · California · USA  
CODAN NORGE AS · Norge  
CODAN TRIPLUS AB · Sverige  
CODAN Limited · Great Britain  
CODAN FRANCE Sarl · France  
CODAN Medical AG · Schweiz  
CODAN ARGUS AG · Schweiz  
CODAN BV · Nederland  
CODAN s.r.l. · Italia  
CODAN Medical GmbH · Österreich  
CODAN Steritex ApS · Danmark  
CODAN Medical ApS · Danmark  
CODAN DEHA ApS · Danmark

### CODAN Product range

- Infusion sets
- Transfusion sets
- Extension lines and manifold connectors
- Infusion and transfusion accessories
- Infusion filters and filter systems
- Neonatology/Paediatric products
- Withdrawal, preparation and administration systems
- CODAN CYTO®
- Chemoprotect® products
- CODAN ENTERAL
- Single use syringes
- Invasive blood pressure monitoring systems
- Infusion pumps
- Continence care products
- Other CODAN Products

Compliance of the established quality management systems with the provisions of EN ISO 13485 and the Council Directive 93/42/EEC has been certified by the relevant, competent notified bodies. Our products are thus entitled to bear the following marks:

### CE 0123 TÜV SÜD Product Service GmbH

CODAN Medizinische Geräte GmbH & Co KG  
23738 Lensahn, Germany

CODAN pvb Critical Care GmbH  
85661 Forstinning, Germany

CODAN US Corporation · Santa Ana, CA 92704, USA  
CODAN ARGUS AG · 6340 Baar, Switzerland

### CE 0543 Presafe Denmark A/S

CODAN Steritex ApS · 3060 Esbjerg, Denmark  
CODAN Medical ApS · 4970 Rødby, Denmark



### Manufacturer

CODAN Medizinische Geräte GmbH & Co KG  
Stig Husted-Andersen Straße 11 · 23738 Lensahn, Germany  
Tel. +49 (0) 43 63 · 51 11 · Fax. +49 (0) 43 63 · 51 12 14  
www.codan.de

The decisive connection

